


Datex-Ohmeda 5330

Reconditioned Agent Monitor

Features

- **Monitors 5 Agents: Enflurane, Halothane, Isoflurane, Sevoflurane, Desflurane**
- **Easy-to-read digital display screen**



Overview of the 5330

These Reconditioned Datex-Ohmeda 5330's are suitable for performing anesthetic vaporizer concentrations checks in the field. These units have been serviced by the OEM and include 5-agent software.



**BC Group International, Inc.
3091 Elm Point Industrial Dr.
St. Charles, MO 63301 USA**

Phone: 314-638-3800 Toll Free: 1-888-223-6763

Fax: 314-638-3200

**Email: sales@bcgroupintl.com
Website: www.bcgroupintl.com**

Datex-Ohmeda 5330

Specifications

General

Height: 4 in, 10.16 cm

Width: 10 in, 25.4 cm

Depth: 16 in, 40.64 cm

Weight: 13.25 lbs, 6 kg

Power Input: 110, 120 VAC, 50/60 Hz

Operational Environment

Temperature: 15 to 35°C

Humidity: 0 to 95% RH, non-condensing

Storage Environment

Temperature: -5 to 50°C

Humidity: 0 to 95% RH, non-condensing

Display

Displayed Parameter: Agent inspired & expired concentration levels (%)

Trendgraph: Last 120 breaths, 20 min, 60 min, 120 min

Sample Flow Rate: 250 mL/min

Response Time: 600 msec

Warm-Up Time: 5 minutes to operation

Compensation: Temperature & barometric pressure

Analog Output: Anesthetic agent waveform

Digital Output: RS-232 compatible

Enflurane, Halothane, Isoflurane, Sevoflurane

Range: 0 to 9%

Accuracy: 0 to 1.0% (+/- 0.2%)

1.1 to 6.0% (+/- 0.3%)

6.1 to 7.0% (+/- 0.35%)

Resolution: 0.1%

Alarms: High expired 0.2 to 10% agent concentration

High inspired 0.2 to 10% agent concentration

Low expired 0 to 9.8% agent concentration

Low inspired 0 to 9.8% agent concentration

Desflurane

Range: 0 to 24%

Accuracy: 0 to 4.9% (+/- 0.2%)

5.0 to 9.9% (+/- 0.5%)

10.0 to 18.0% (+/- 1.0%)

18.1 to 22.0% (+/- 2.0%)

Resolution: 0.1%

Alarms: High expired 0.2 to 10% agent concentration

High inspired 0.2 to 10% agent concentration

Low expired 0 to 9.8% agent concentration

Low inspired 0 to 9.8% agent concentration

BC Group International, Inc.
3091 Elm Point Industrial Dr.
St. Charles, MO 63301 USA

Phone: 314-638-3800 Toll Free: 1-888-223-6763

Fax: 314-638-3200

Email: sales@bcgroupintl.com

Website: www.bcgroupintl.com

