


# 7245 Sterile Petri Dishes



## Description

### Unique patented design for easy handling and storage

- Easy to use. Opens easily with one hand, yet close to a tight seal.
- Use less space on the lab bench or in the incubator with easy stacking base.
- Gamma irradiated for sterilization. No dangerous EtO residuals.
- Available with or without absorbent pads.
- Bulk packaging for additional value.

## Applications

- Ideal for microbiological analysis when performing the Membrane Filter (MF) Technique.
- Petri dishes with absorbent pads can be used with broth media or pour agar into dishes without absorbent pads.



## Specifications

### Petri Dish

#### Material of Construction

Polystyrene

#### Dimensions

Height: 9.0 mm (0.35 in.)

Diameter: 50.0 mm (1.97 in.)

#### Filter Size

Accepts 47 mm membrane filter

#### Sterilization

Gamma irradiation

### Absorbent Pads

#### Composition

Cellulose

#### Typical Thickness

0.9 mm (35 mils)

BC Group International, Inc.  
3081 Elm Point Industrial Dr  
St. Charles, MO 63301 USA  
314-638-3800  
314-638-3200 Fax  
1-888-223-6763



[www.bcgrouptl.com](http://www.bcgrouptl.com)